



5
in 1

multifunctional
Oxygen
Facial Machine



·Improve acne

·Increase flexibility

·Deep moisturizing

·Oxygen revitalizing

·whitening

·wrinkle removal

·skin lifting




Oxygen facial machine is designed to revitalize the skin. The professional oxygen facial machine combines oxygen spray gun, Bipolar RF, Hydradermabrasion, diamond dermabrasion and PDT mask. It complements all skin care products and facial services and can be easily incorporated into your existing beauty salon treatment menu. Age, pollution and sun exposure all contribute to the degeneration of skin capillaries and the portable oxygen facial machine delivers a pressurized flow of concentrated oxygen and active cosmetic serums to re-oxygenate and revitalize the skin. Lightweight, portable, low maintenance and whisper quiet, it can be an excellent addition to your treatment offerings and because it's comfortable.

This portable oxygen facial machine PL-6 creates instant beautiful healthy skin results. **Clients and/or patients obtain instant gratification after their first treatment.**

Working Principle

Nano-oxygen lance (deep hydration) nanoparticles spray oxygen, water particles are fine, directly penetrate the muscles bottom, 16 times the ordinary spray effect.




Oxygen spray gun



Hydra dermabrasion

Hydradermabrasion utilizes the natural healing powers of water and oxygen to effortlessly exfoliate the skin without the use of hard crystals or abrasive textured wands, generating deeply hydrated and healthier looking skin.


Deep cleansing·
vacuum absorption·
fundamentally clean without damaging the skin·



The extraction of oxygen from the air which is compressed into an active factor by a molecular filter, and the release of oxygen into the skin, is an effective and safe method of delivering oxygen to the skin to remove wrinkles and smooth skin.



Diamond Peel




Radio frequency

The radio waves enter the deep part of the skin, and the internal natural resistance movement causes the skin to heat up. The use of dermal collagen cellulose immediately produces a contraction principle at 45~60°C, which can smooth wrinkles, tighten the skin and delay skin aging.

PDT wavelength

PDT emits light of a specific wave length, penetrates the surface of the skin, reaches the subcutaneous fat layer, stimulates capillary blood flow, increases cell production, increases blood circulation, and produces rich collagen and elastin. Achieve easy acne, remove wrinkles, increase skin elasticity and restore skin vitality.



Seven color spectrometer
PDT mask
(Seven spectral free switching)



Red light repair skin



White light shrinks pores



Yellow light fade melanin



Blue light hyperplasia collagen, firming skin



Green light balance oil



Purple light whitening skin

Technology

Oxygen spray gun probe

Put the nutrient solution into the work head, let the product spray directly on the skin, let the skin absorb better, inject oxygen and hydrate, activate cells.

Supply nutrients and moisture for the skin rapidly, smooth skin, lightening stains, improve chlorosis, dull conditions, make skin rosy, white, transparent, tender, tightening indeed, with special spa oils will be better.



Hydra facial water dermabrasion

Hydra facial- water dermabrasion is an exciting NEW technology that combines the efficacy of Microdermabrasion, vacuum system and new Aqua Fuse hydration system. It is a lot gentler than crystal microdermabrasion or diamond dermabrasion as only water is used, and we use professional quality oxygen facial machine.

Hydra-dermabrasion (or hydra microdermabrasion) applies mechanical and chemical peeling simultaneously. A hydra-dermabrasion machine includes a water jet, an air compressor, a two-way flow control valve, a purified container and a waste water tank.



Diamond dermabrasion probe

The diamond microdermabrasion (crystal-free) system replaces the micro-crystal flow with a diamond tip wand. Areas closer to the eyes and mouth can be abraded more accurately with diamond abrasion – no risk of stray crystals. Provides a Non-Surgical skin resurfacing procedure by using sterile diamond heads to peel and rub off the dead cells at the top skin layer followed by vacuuming/suction removing particles along with any dirt and dead skin.

This procedure removes skin debris, scars, blemishes, wrinkles and uneven pigmentation of the skin. Moderate usage of diamond dermabrasion accompanied with skin products that penetrate through the top skin layer and reach the dermis level help to replenish natural nutrients and restore cell activity promoting healthy skin.



BIO Handle

The BIO working handle has calming, smoothing and firming effect, shrink pores, and reduce the swelling excrete by the process of metabolism, It is has great effect on skin lifting and the vessel and nerve will not be damaged.



PDT mask

Praised as “Magic light”, Photon dynamical instrument is a biological active cold light, and not produce high heat. It applies low-energy and unique optical technology to skin care and treatment, natural and gentle; has prominent effects, low fatalness and no side effect. It's skin nurse specialist, is fit for all kinds of skin, especially useful for seriously irritated skin, acne skin and has notable effect for those in sub healthy state to recover and do skin protection and treatment.

Photon dynamical nurse therapy is a new technology, using the photon power to do skin maintenance and treatment, transfer light energy to cell energy, then accelerate cell growth and blood circulation, stimulate fiber cells to produce collagen proteins, increase skin elasticity and take important part in rehabbing aging and acne skin, lightening flecks, skin upgrade and tighten and easing burnable skin.



Technical Parameters

Max output	250VA
Voltage & frequency	VAC 110-220;50/60Hz
Screen	7" Touch screen
Power	250 W
Package	47*62*63cm
Vacuum	90 KPa
Technology	Hydro dermabrasion, PDT Mask
Working handles	Hydra dermabrasion, Diamond dermabrasion, BIO, Spray gun. PDT Mask
Net weight	20/35KG

System Benefits

- 7-inch Touch screen.
- Updating 5 types of operational probes.
- Professional crystal heads.
- Hydro-dermabrasion, applicable to regular or sensitive skin, or skin with whitehead, comedo, acne, etc.
- Cleaning & washing – deep cleaning, clear skin status, remove blackhead, remove deep skin dirt.
- Effective and direct moistening – supply sufficient water molecules to skin while cleaning.





Whatsapp: 0086 15001300057

Web: <https://www.prettylasers.com>

Email: info@prettylasers.com

Phone: +0086 15001300057



Contact With
PrettyLasers