



## Body Sculpting Machine

- Coolsculpting
- Vacuum
- RF
- Shock wave
- Cavitation
- Lipolaser

Two cryo handles work at the same time.

Coolsculpting machine for sale is an advanced equipment using the newest coolsculpting technology (the same as cryolipolysis) to noticeably eliminate fat without surgery. Cool sculpting is a non-invasive procedure, which uses a targeted cooling process to kill the fat cells underneath the skin, freezing them to the point of elimination.

Only fat cells are frozen. Your healthy skin cells remain, well, healthy. No knives. Non suction hoses. No needles. No scars. Once crystallized, the fat cells die and are naturally eliminated from your body. Especially the belly fat, the effect is very obvious.

The cool sculpting machine is indicated for cold-assisted lipolysis (breakdown of fat) of the submental area, thigh, abdomen, and flank, or "love handles" in individuals with a Body Mass Index (BMI) of 30 or less. Cool sculpting is the latest new technology for body slimming, body shaping, body contouring, weight loss, belly fat reduction, cellulite reduction etc.

CoolSculpting uses a procedure known as cryolipolysis. The cryo handle can absorb the fat in the treatment area to the cooling plate by a vacuum pressure, which cools the fat to a freezing temperature. And the treatment probe can control the cooling temperature precisely during the treatment, which can eliminate fat cells in a particular part of the body, such as the belly fat. When the fat cells are in the cooling environment which be precisely controlled, it will begin to do natural decomposition and clearance, to make the fat layer become thinner gradually.

# Product Principle




**Cryo Handle**



CoolSculpting uses a procedure known as cryolipolysis. The cryo handle can absorb the fat in the treatment area to the cooling plate by a vacuum pressure, which cools the fat to a freezing temperature. And the treatment probe can control the cooling temperature precisely during the treatment, which can eliminate fat cells in a particular part of the body, such as the belly fat. When the fat cells are in the cooling environment which be precisely controlled, it will begin to do natural decomposition and clearance, to make the fat layer become thinner



**LipoLasers**




Laser for fat reduction, use the 660nm and 980nm soft laser to destroy the deep fat cells, broke the cell membrane, change the fat to tri-glyceride, then leave out the body, thanks to metabolism.



**Shock Wave**




The shock wave at a frequency of 1-16HZ up to 8cm under the skin, causing strong resonance rotation of each tissue cell, local heat increase, and the instantaneous temperature rise reaches 43 °C. Promotes blood flow, stimulates damaged cells, activates collagen cells, and promotes the reorganization of fibrous connective tissue.




Using strong Ultrasound to shake fat cells, reducing fatty acids, enhancing blood circulation, and finally letting fatty acids excrete through lymphatic drainage, thereby reducing body weight. It can also lighten cellulite and firm skin.

**Cavitation**



Radio frequency radio waves enter the deep part of the skin, and the internal natural resistance movement causes the skin to heat up. The use of dermal collagen cellulose immediately produces a contraction principle at 45-60C', which can smooth wrinkles, tighten the skin and delay skin aging.

**Three polar RF**



**Vacuum RF**

A vacuum elevates the target tissue, bringing it closer to the energy source. Directly heats the targeted areas, Bi-polar radio frequency (RF) energy safety guides heat throughout the targeted area, so that it heals faster without damaging the skin.

## CoolSculpting Machine for Sale Treatment Handles

**LONGER COOLING PLATES** – The maximum treatment area is over 40% of the cooling plate.

**MEDICAL SILICA GEL** – Effectively observe the treatment and protect the skin.

**TOUCH SCREEN** – Real-time monitoring treatment head temperature change error alarm.

**UNIQUE TREATMENT HEAD** – Unique double chin treatment handpiece.

**LOWEST TEMPERATURE  
CAN REACH -16°C**  
full ice freeze on the entire cooling plates



°C  
- 16.7

Actual test  
temperature

This professional cryolipolysis machine features with 5 handles for body sculpturing, five treatment is by fat freezing application, including 360degree cryo handle, two small handles, mid handle, big

### The Largest Applicator

This is the largest applicator used to treat belly hangover and larger areas on the abdomen, it allows for maximum fat reduction. This applicator can be used on the arms, hips, and men's chest area.



**The Largest Applicator: treat area (250mm \*120mm) for belly, back, buttock**

---

### The Medium Applicator

Medium Applicator is great for those annoying flanks. So, this applicator on how it curves around your sides to adds a more defined waist.



**The Medium Handle: treat area (200mm \* 100mm) – for waist, thigh,**

---

### The Small Applicator

This is excellent for underarm fat blasting. This applicator is great for those tiny pockets of fat that won't go away around the knees, arms, and back of thighs. Sculpting banana rolls and bra bulges has never been so simple.



**Small Handle: treat area (120mm \* 60mm) – for thigh, arm, crus**

---

### 360° Mini Cryo Handle

The 360° Mini Cryo Handle is the smallest applicator. we're getting amazing results for double chin fat reduction with this applicator. Combine with three polar RF for non-surgical neck lift and fat reduction. Three polar RF tightens the neck, while the 360° Mini Cryo Handle removes the fat.



**360° Mini Cryo Handle: treat area (76mm \* 60mm) – for double chin,**



### Cavitation Handle

The cavitation handle also plays a role in shaping your target areas, “chiseling” and breaking down tissue to give you a slimmer appearance.

Remove fat  
Lose weight  
Full body available.

---

### Three Polar RF Handle

RF handle can deep direct fat body, has the directional positioning RF output excellence, make fat cells organizations in fast active state, produce the cell body heat of the friction heat, in order to achieve the effect of dissolved fat.

Body slimming  
Lifting contours  
Tightening skin



---

### Shock Wave Probe

Shockwave is delivered by placing a probe up against the areas of cellulite, as it passes over the area it emits Radial Waves through the skin.

Meridian massage  
Pain point physiotherapy  
Deep penetration



---

### LipoLaser

The laser lipolysis opens and releases glycerol, water and free fatty acids, thus, causing to lose weight and inches.

Reduce cellulite  
Relieve physical fatigue



---

### Vaccum RF Handle

The Vacuum RF handle uses suction and pressure to release the fluids in the fat cells that cause cellulite.

Dredge the meridians  
Speed up metabolism  
Promote absorption



NEW DESIGN NEW LOOK

# FUNCTIONS

## SKIN TIGHTEN

Lifts and tightens skin  
restores skin elasticity and  
rejuvenates skin

## LOSE WEIGHT FAST

Slimming and shaping  
Painless, non-invasive, no  
side effects, obvious effect

## STIMULATE LYMPHATIC DETOX

Promote metabolism  
activate subcutaneous blood circulation

## REDUCE CELLULITE

Smooth skin  
Tighten skin, relieve skin  
sagging, reduce orange peel

## PHYSICAL RELAXATION

Relieve physical fatigue  
Meridian massage, improve  
blood circulation, stimulate  
acupuncture points





# PRODUCT

## >>> DETAILS



Vacuum RF



Cavitation



Three polar RF

Screen



Lipolaser

Switch button




Cryo handle

(Two cryo handle work at the same time )



Shock wave handle



Diffuser

Whatsapp: 0086 15001300057

Contact With  
PrettyLasers

Web: <https://www.prettylasers.com>

Email: [info@prettylasers.com](mailto:info@prettylasers.com)

Phone: +0086 15001300057

