

**5Handles
Professional
Cryolipolysis
Slimming
Machine**

PL-CT80



5 Handles professional cryolipolysis slimming machine is a non-invasive procedure, which uses a targeted cooling process to kill the fat cells underneath the skin, freezing them to the point of elimination. This professional machine features with 5 handles for body sculpting, five treatment is by fat freezing application, including 360degree cryo handle, two small handles, mid handle, big handle.

Four cryo handles can work at the same time.

Professional electronic structure to ensure stable working.

CoolSculpting Machine for Sale Treatment Handles

This professional cryolipolysis machine features with 5 handles for body sculpting, five treatment is by fat freezing application, including 360degree cryo handle, two small handles, mid handle, big handle.

LONGER COOLING PLATES – The maximum treatment area is over 40% of the cooling plate.

MEDICAL SILICA GEL – Effectively observe the treatment and protect the skin.

TOUCH SCREEN – Real-time monitoring treatment head temperature change error alarm.

UNIQUE TREATMENT HEAD – Unique double chin treatment handpiece.



Anti-freeze membrane

Anti-freeze membrane is an important and necessary part of every, single treatment with the most advanced method of eliminating fat - Cryolipolysis, which would provide ideal guidance cold with minimum loss of energy and also acted as a better protection of the surface layer of the skin.

Effectively protect the skin from direct exposure to the low temperature.

Enables rapid management of cold air into the subcutaneous tissue to stored fat cells in order to affect the most efficient Cryolipolysis.

MANY SIZES OF ANTIFREEZE MEMBRANE



SIZE:12*12	WEIGHT:30g	1 pc/bag
SIZE:28*28	WEIGHT:100g	2 pc/bag
SIZE:32*32	WEIGHT:110g	2 pc/bag
SIZE:34*42	WEIGHT:70g	1 pc/bag

SIZES FOR FREEZING HEADS

Four Cryo handles can work at the same time

LONGER COOLING PLATES

The maximum treatment area is over 40% of the cooling plate



MEDCAL SILICA GEL

Effectively observe the treatment and protect the skin



TOUCH SCREEN

Real-time monitoring
treatment head temperature
change error alarm



UNIQUE TREATMENT HEAD


Unique double chin
treatment handpiece



Big handle: treat area
(250mm*120mm)
—for belly, back,
buttock etc.



Mid handle: treat area
(200mm*100mm)
—for waist, thigh, etc



Small handle treat area
(120mm*60mm)
—for thigh, arm, crus



360° Mini cryo handle
treat area
(76mm*40mm)
—for double chin, arm



Technical Parameters

Voltage	220V/50Hz or 110V/50Hz or 110V/60Hz
Cryo output temperature	5 degree~ -16°C degree
Vacuum output pressure	0-100Kpa
Cooling liquid	Pure water, with two 9L
Default working time	0-60 minutes
Power consumption	≤400W
Weight	45KGS (NW) / 65KGS (GW)
Machine Measure	105cm*470*550mm
Outside Package	59*57*131mm

Advantages

Power supply, imported from Meanwell, Taiwan, to ensure the machine to keep working over 12 hours.

Electromagenatic valve, imported from USA, to control the suction and release precisely and stably.

Cooling system Water cooling + Wind + Semiconductor.

Contact Skin Parts: Medical use silicon, soft to make customer feel comfortable during treatment.

Proportional valve, imported from Germany, to control the vacuum power.

Four Air pump, imported from Taiwan, Four drainage system, double water pump.

Cooling temperature: -16~5 degree.

Vacuum Power: 100KPa(max).

Quality ABS case

Four cooling system with 14L water tanks which ensure machine can working non-stop 12 hours etc.

Four Cro handles can work at the same time. Professional inner electronic structure to ensure stable working.

THREE COLOR

WHITE



GREEN



BLUE



Whatsapp: 0086 15001300057

Web: <https://www.prettylasers.com>

Email: info@prettylasers.com

Phone: +0086 15001300057

Contact With
PrettyLasers

